Design and Fabrication of Micro-Tip Arrays for Nano-Scale Atom-Probe Tomography

Robin Peter, University of Chicago NNCI Site: Northwestern University Project Mentors: Dr. Ying Jia, Dr. Dieter Isheim Principal Investigator: Dr. David Seidman







Northwestern

Universitv

Y Overview: Atom-Probe Tomography (APT)¹



Applications:



Image sources: from left to right, see References [2], [3], [4].



UNIVERSITY OF



Schematic of APT Mechanism (Seidman⁵)

- Laser/voltage pulse generates electric field
- Field evaporation of atoms as ions detected, measured, analyzed
- Requires needleshaped specimen

Specimen preparation is the major limitation for wide-spread APT use.







Northwestern

Y

Specimen Preparation Factors

Functional Factors	Parameters ⁶	
Electric field strength	30-50 V nm⁻¹	
Needle tip radius	20-150 nm	
Semi-angle taper	< 10 [°]	
Material resistivity	< 0.05 Ω cm	

Other Factors	Options		
Throughput	Single needles	Micro-arrays	
Shape	Pre-sharpened micro-tips (PSMs)	Flat-top micro-posts (FTM) ⁷	

CHICAGO



Northwestern

Y Reconsidering Commercial Arrays



Broken micro-tip arrays.

Information source: communication with Dieter Isheim, Northwestern University (July 2017).

- **1. Cost:** \$250 / 22-tip array
- 2. Design: marker posts for tip identification interfere with specimen visibility
- **3.** Fragility: tips are unprotected in the event of manual mishandling of an array
- 4. Flexibility: customization beyond commercial manufacture may allow joint analysis between APT, TEM, and other characterization equipment







Northwestern

PROJECT OBJECTIVE: Fabricate a microarray satisfying APT parameters.



SEM images of posts cut with (a) reactive-ion etching and (b) a dicing saw, from Kelly et al.³







Northwestern University

Design: Two Tiers

Total Height: 150-180 μm **Top Layer:** FTM/needles

- 40 µm height
- ~2 μm tip platform
 Base Layer: Angular structure
- 100-200 μm sides
- 110-140 μm height



Process Overview:

- 1. Fabricate needles
- 2. Fabricate base
- 3. Ensure needles are protected during (1.) and (2.)







HYNE Soft and Hybrid Nanotechnology Experimental Resource





METHODS AND RESULTS

SHYNE Soft and Hybrid Nanotechnology Experimental Resource









Fabricating: Needles

Multiple etching processes generate the desired needle geometry.







Northwestern

Method 1: Wet KOH Etching

Goal: sharpen cylinders into FTMs or needles.





SEM image of KOH-etched Si post, 4448X.







Method 2: Isotropic DRIE

Goal: sharpen cylinders into FTMs or needles.





SEM image of DRIE on 15- μ m Si post, 4695X.







THE UNIVERSITY OF CHICAGO





Method 3: Reverse-Order DRIE

Line Profiles: Effect of Multi-step DRIE on Silicon Wafer



'HE UNIVERSITY OF

CAG

Soft and Hybrid Nanotechnology SH Experimental Resource

Needle Structures: Conclusions

- Wet chemical KOH etching is not viable for producing needle-like structures at the desired scale
- Isotropic DRIE is viable for producing a needleshaped structure, but will not result in the desired height without further processing

Looking ahead:

 Continue to explore deep-trench and isotropic DRIE parameters

5HYNE Soft and Hybrid Nanotechnology Experimental Resource





Northwestern

Jniversity

Fabricating: Base Structures

Simulation: KOH Etch of Base Structures



Fig. 9: Simulation of 30% KOH wet etch on Si using Anisotropic Crystalline Etching Simulation (ACES) software.⁹

Parameters: <100> face and <110> edge wafer 5-minute increments at 1.4 µm/min

(Note: tested etching rate for the procedure of this project is approximately half the simulation rate.)

SHYNE Soft and Hybrid Nanotechnology Experimental Resource





Northwestern

Fabricating: Base Structures

- Etching rate of 30%
 KOH soln. at 85°C:
 ~40 μm/hr
- 3 hrs of KOH etching will produce a 120-µm structure

Line Profiles: Effect of KOH Etching on 100-µm Square-Array Mask



→ No Etch → 30 min → 60 min

HE UNIVERSITY OF

Northwestern



Fabricating: Base Structures

3 Hr KOH Etch (100 μm square) Over-etched

SH

2 Hr KOH Etch (100 μm square) Correct geometry





THE UNIVERSITY OF CHICAGO





Base Structures: Conclusions

- Wet chemical KOH etching is viable for producing large structures to support needles
- **3 hours of etching** (under stated conditions) will produce the desired height

Looking ahead:

• Use larger mask size to retain platform









Test: Needle Protection

Q: Is SiN masking effective at protecting upper needle structures during KOH etching of the base?



Deposit SiN



Add, spin-coat photoresist



Maskless aligner exposure



Rest Postbake



Develop photoresist



RIE











Test: Needle Protection

Result: adequate protection, inadequate photoresist











Northwestern

Conclusions and Further Research

Microfabrication techniques can produce components of a **multi-tiered APT micro-array**



CHICAGO



Northwestern



Acknowledgements

This work utilized Northwestern University Micro/Nano Fabrication Facility (NUFAB), which is partially supported by Soft and Hybrid Nanotechnology Experimental (SHyNE) Resource (NSF ECCS-1542205), the Materials Research Science and Engineering Center (NSF DMR-1121262), the State of Illinois, and Northwestern University.

This work made use of the EPIC facility of Northwestern University's NUANCE Center, which has received support from the Soft and Hybrid Nanotechnology Experimental (SHyNE) Resource (NSF ECCS-1542205); the MRSEC program (NSF DMR-1121262) at the Materials Research Center; the International Institute for Nanotechnology (IIN); the Keck Foundation; and the State of Illinois, through the IIN.

The local-electrode atom-probe tomograph at the Northwestern University Center for Atom-Probe Tomography (NUCAPT) was acquired and upgraded with equipment grants from the MRI program of the National Science Foundation (NSF DMR-0420532) and the DURIP program of the Office of Naval Research (No0014-0400798, No0014-0610539, No0014-0910781). NUCAPT received support from the MRSEC program (NSF DMR-1121262) at the Materials Research Center, the SHyNE Resource (NSF NNCI-1542205), and the Initiative for Sustainability and Energy at Northwestern (ISEN).

Soft and Hybrid Nanotechnology Experimental Resource











References and Image Sources

- 1. Kelly TF, Larson DJ, Thompson K, Alvis RL, Bunton JH, Olson JD, Gorman BP. 2007. Atom Probe Tomography of Electronic Materials. *Annu. Rev. Mater. Res.* 37:681-727.
- 2. Crusito, RL. n.d.. Microchip integrated into motherboard [online image]. Retrieved from www.azonano.com/article.aspx?ArticleID=3978.
- 3. AMG Advanced Metallurgical Group. 2012. AMG Superalloys slideshow photo [digital photograph]. Retrieved from www.amg-s.com.
- 4. CVD Equipment Corporation. 2017. Silicon nanowires SEM 10 K X [digital photograph]. Retrieved from www.firstnano.com/portfolio-item/silicon-nanowires-sem-10-k-x/.
- 5. Seidman, DN. 2007. Three-dimensional atom-probe microscopy: advances and applications. *Annu. Rev. Mater. Res.* 37:127-158.
- 6. Kelly TF, Miller MK. 2007. Atom probe tomography. *Review of Scientific Instruments* 78:031101.
- 7. Nayak AP, Logeeswaran VJ, Islam MS. (n.d.). *Wet and Dry Etching* [PDF file]. Retrieved from web.ece.ucdavis.edu/~anayakpr/Papers/Wet%20and%20Dry%20Etching_submitted.pdf
- 8. Thompson K, Lawrence DJ, Larson DJ, Olson JD, Kelly TF, Gorman BP. 2007. In-situ site-specific specimen preparation for atom probe tomography. *Ultramicroscopy* 107:131–39.
- 9. Zhu Z, Liu C. 1998. Anisotropic crystalline etching simulation using a continuous cellular automata algorithm.





Northwestern